

high performance plastics


# STAINLESS STEEL ACCESSORIES

Perfectly fitted — for shining performances

The PurQuip® stainless steel accessories are designed to fulfil highest requirements and impress through excellent materials and workmanship. High performance plastic and polished stainless steel surfaces guarantee GMP and cleanroom conformity along with long-term resistance to steam sterilization (autoclaving).

- ▶ Mirror polished stainless steel surfaces
- ▶ Best processing quality and conformity
- ▶ Durable, autoclavable
- ▶ Ergonomic quick-change system
- ▶ All frames and handles can be combined with each other

electropolished stainless steel surface



## PurQuip® KEU40 Folding mop frame

40 cm flat mop system

- ▶ Stainless steel, electropolished
- ▶ Handle socket made of special plastic
- ▶ Autoclavable, disinfectant resistant
- ▶ **Quick-change system with clip**
- ▶ **Combinable with all PurQuip® handles**
- ▶ **Allows mop pick-up without hand contact**

Art.-No. 2151003

A ISO GMP 3 4 5 A/B 6 7 C 8 D



## INNOVATION

### PurQuip® quick-change system

Allows for the quick, simple exchange of the handle and mop frame without tools. The special clip system is not only very diverse, but also simple: All PurQuip® mop frames and handles can be connected with each other in one action.



### Touchless mop system

In combination with the stainless steel mop frame PurQuip® KEU40, all PurMop® flat mops can be attached, wetted, wrung out and detached without hand contact. The touchfree concept guarantees safe, hygienic cleaning processes without cross-contaminations, and with the best possible ergonomics.



# MOP SYSTEMS

Stainless steel accessories

**TIP**  
All PurQuip® frames and handles are compatible!



## PurQuip® TE30 PurQuip® TE40

**2 widths**  
30 cm | 40 cm

Mop frame, T-Form

- ▶ Stainless steel, electropolished
- ▶ Autoclavable, disinfectant resistant
- ▶ Light-weight, simple handling
- ▶ Quick-change system with clip
- ▶ Widely adjustable (180° pivoting)
- ▶ For the cleaning of walls, ceilings and equipment
- ▶ Suitable for mops R32 | R42 (page 21)

TE30: 30cm | TE40: 40cm  
Art.-No. TE30: 2161001 | TE40: 2161002

**A** ISO GMP 3 4 5 A/B 6 7 C 8 D



PurQuip

## PurQuip® SEP50 PurQuip® SEP140

**2 lengths**  
50 cm | 140 cm

Handle, rigid

- ▶ Stainless steel, electropolished
- ▶ Autoclavable, disinfectant resistant
- ▶ Plastic handle for easy handling
- ▶ Quick-change system with clip
- ▶ Reduced weight

SEP50: 50cm | SEP140: 140cm  
Art.-No. SEP50: 2131013 | SEP140: 2131003

**A** ISO GMP 3 4 5 A/B 6 7 C 8 D



PurQuip

## PurQuip® SET2070 PurQuip® SET2100

**2 lengths**  
2 x 70 cm | 2 x 100 cm  
individually adjustable

Telescopic handle, length adjustable (2-piece)

- ▶ Stainless steel, electropolished
- ▶ Autoclavable, disinfectant resistant
- ▶ Plastic handle for easy handling
- ▶ Quick-change system with clip
- ▶ Can be adjusted individually to body size
- ▶ Reduced weight
- ▶ Specially suitable for wall and ceiling cleaning (overhead work)

SET2070: 2 x 70cm (140cm)  
SET2100: 2 x 100cm (200cm)  
Art.-No. SET2070: 2131006 | SET2100: 2131005

**A** ISO GMP 3 4 5 A/B 6 7 C 8 D



PurQuip



MOP ACCESSORIES

# ISOLATOR CLEANING TOOL

Mop accessories, stainless steel



**Set:**  
SEP50 (2131013) and  
MEP20 (2151006)

## PurQuip® ICT2050 Isolator cleaning system

20cm ICT system

- ▶ Compact and handy cleaning tool for disposable mop heads
- ▶ Stainless steel, electropolished
- ▶ Autoclavable, disinfectant resistant
- ▶ GMP-compliant materials for critical application
- ▶ **Quick-change system with clip**
- ▶ **Adjustable mop frame**
- ▶ **Compatible with PurMop® EC20-S / EC20-SV**
- ▶ For cleanroom isolators, workbenches and equipment

Art.-No. 2131015

**A** ISO GMP 3 4 5 A/B 6 7 C 8 D



Optional:  
presaturated, ready for use

## PurMop® EC20-S | [EC20-SV] Cleanroom disposable mop, sterile, [sterile presaturated]

20 cm ICT system | 20x 2 pieces per case

- ▶ 100% polyester yarn
- ▶ Knitted, embossed material, multi-layer
- ▶ Extra fine wetting properties and low particle emission
- ▶ Particularly suited for stainless steel and glass surfaces
- ▶ **Pre-soaked with sterile isopropanol / WFI water 70/30 (EC20-SV)**
- ▶ **Pocket mop for firm fit on the frame**
- ▶ **Sterile, double packed**
- ▶ For cleanroom isolators, workbenches and equipment

Art.-No. EC20-S 2111111 | [EC20-SV 2111112]

**S** ISO GMP 3 4 5 A/B 6 7 C 8 D



## INNOVATION

### The PurQuip® Isolator cleaning system

The new cleaning tool was developed for simple and reliable cleaning of isolators, RABS systems, workbenches and equipment in the cleanroom. The special version with a short handle and pivoting mop frame with adjustable fixture was designed with high user friendliness in mind.

### Maximum safety

With the use of the ICT2050 the risk of contamination is significantly reduced – cleaning of isolators is conveniently performed from the outside without danger for the user. The dual articulated joint of the mop frame allows access everywhere. The special disposable mops are made of multi-layer polyester material to ensure maximum cleanroom suitability – crucial for cleaning of critical surfaces! Two validated models are available: dry or presaturated with IPA 70/30. Both are delivered sterile, in double packaging and ready to use.



length adjustment by the push of a button

robust plastic

ergonomic handle

# ALUMINIUM/PLASTIC ACCESSORIES

## The ease of cleaning

The PurQuip® plastic accessories impress through clever technics, light-weight and the best ergonomics along with a large selection of handles in various lengths. Special smooth plastic surfaces facilitate an easy cleaning and good resistance.

- ▶ Ergonomic and light-weight
- ▶ Large selection of handles in different lengths
- ▶ Clever adjustment and innovations
- ▶ Smooth plastic and aluminium surfaces

### INNOVATION

#### PurQuip® SAK170 Telescopic handle

- ▶ Length adjustment by the push of a button (single hand operation)
- ▶ Can be adjusted individually to the cleaning task/user
- ▶ Grip-ball for optimal ergonomics and secure hold



Optimal for walls and ceilings!


#### PurQuip® KKV40 Folding mop frame

- 40 cm System
- ▶ Polypropylene, smooth and resistant
  - ▶ Hands-free opening mechanism
  - ▶ Autoclavable, disinfectant resistant
  - ▶ Quick-change system with clip
  - ▶ Locking device for simple ceiling cleaning
  - ▶ Light-weight and ergonomic

Art.-No. 2151005

A ISO GMP 3 4 5 A/B 6 7 C 8 D

PurQuip



#### PurQuip® SAK170 Telescopic handle, length adjustable (2-piece)

- 170 cm length
- ▶ Aluminum main body encased in plastic
  - ▶ Single hand adjustment
  - ▶ Quick-change system with clip
  - ▶ Grip-ball for comfortable handling
  - ▶ Light-weight and robust

Art.-No. 2131007

A ISO GMP 3 4 5 A/B 6 7 C 8 D

PurQuip





# MOPP SYSTEMS

Aluminium / plastic accessories



## PurQuip® SAS140 A

Handle, rigid

140 cm length

- ▶ Aluminium main body encased in plastic
- ▶ Autoclavable, disinfectant resistant
- ▶ Grip-ball for comfortable handling

Art.-No. 2131010

ISO 3 4 5 6 7 8  
GMP A/B C D

PurQuip



## PurQuip® SAK80

Telescopic handle, length adjustable (2-piece)

80 cm length

- ▶ Aluminium main body encased in plastic
- ▶ Disinfectant resistant
- ▶ Grip-ball for comfortable handling
- ▶ Ultra short, specially for cleaning of equipment and workbenches

Art.-No. 2131008

ISO 3 4 5 6 7 8  
GMP A/B C D

PurQuip



## PurQuip® SA180

Telescopic handle, length adjustable (2-piece)

180 cm length

- ▶ Aluminium
- ▶ Disinfectant resistant
- ▶ Quick-change system with clip

Art.-No. 2131009

ISO 3 4 5 6 7 8  
GMP A/B C D

PurQuip



## PurQuip® SA400

Telescopic handle, length adjustable (3-piece)

400 cm length

- ▶ Aluminium
- ▶ Disinfectant resistant
- ▶ Quick-change system with clip
- ▶ Specially suitable for overhead and ceiling cleaning

Art.-No. 2131012

ISO 3 4 5 6 7 8  
GMP A/B C D

PurQuip



Extra long!

MOP ACCESSORIES