



Comfortably Protected

The optimal protective products not only offer us uncompromising safety, but also ensure maximum wearer comfort. This allows us to focus on the essential – the manufacturing of our products.

PROTECTIVE EYEWEAR

A clear vision for safety

The protection of the cleanroom operator is just as much a task of protective eyewear, as is the protection of the clean environment. Discover PurGuard™ products that will enhance your vision and keep you focused on the job at hand.

PROTECTIVE EYEWEAR	74
▶ Cleanroom goggles	76

S Sterile **A** Autoclavable



Note:
There are approximately 1 million germs per cm² on the human forehead and scalp!

CLEANROOM GOGGLES

A safe choice

Particularly from the area around the eyes, which is still often left uncovered, countless particles are released every minute from which your product needs to be protected. Only complete protection of the head and face area with special masks and cleanroom suitable goggles offers uncompromising safety in the cleanroom.

- ▶ For work in the cleanroom and aseptic environment
- ▶ Disposable products with comprehensive validation
- ▶ Fulfills the requirements of DIN ISO and GMP
- ▶ For combined product and personnel protection

PurGuard™ GOGGLES

At a glance

REUSABLE


- ▶ Long term reusable
- ▶ Economic and resource saving
- ▶ CE certified for personal protection
- ▶ High airflow through special ventilation
- ▶ Comfortable fit
- ▶ Different models available



DISPOSABLE

- ▶ Delivered sterile
- ▶ Safety through full validation
- ▶ CE certified for personal protection
- ▶ Anti-fog coated lens
- ▶ Can be worn over prescription glasses
- ▶ Clear frame for bright vision

LIFE CYCLE STEAM STERILIZATION



Reusable models

CLEANROOM GOGGLES



All PurGuard™ goggles are CE certified according to DIN EN 166:2002!



Suitable for wearing over prescription glasses

Panorama lens for allround vision

PurGuard™ SV-900 panorama Autoclavable goggles

- 10 pieces per unit (individually packed)
- ▶ Soft, very comfortable frame material
- ▶ Robust lens
- ▶ Suitable for 20 autoclave cycles (at 121°/134° C).
- ▶ Adjustable silicone headband (latex-free)
- ▶ **Panorama shape for allround vision**
- ▶ Rotatable strap hinges allow for broad adjustability
- ▶ Special ventilation concept prevents fogging: direct on the top, indirect on the bottom



Art.-No. 2962003



PurGuard™ SV-800 Autoclavable goggles

- 10 pieces per unit (individually packed)
- ▶ Soft, very comfortable frame material
- ▶ Robust lens
- ▶ Suitable for 20 autoclave cycles (at 121°/134° C).
- ▶ Adjustable silicone headband (latex-free)
- ▶ **Specially suitable to be worn over prescription glasses**
- ▶ Spacious, extra wide design
- ▶ Special ventilation concept prevents fogging: direct on the top, indirect on the bottom



Art.-No. 2962004



PurGuard™ EV-100 Sterile disposable goggles

- 50 pieces per unit (5 bags of 10 pieces, individually packed)
- ▶ Ultra soft, comfortable frame material
- ▶ Robust lens with special anti-fog and anti-scratch coating
- ▶ Double packed, sterile
- ▶ Adjustable silicone headband (latex-free)
- ▶ **Clear frame material and extra wide design for prescription glasses**
- ▶ Special ventilation concept prevents fogging: direct on the top, indirect on the bottom
- ▶ No sterilization required after usage



Art.-No. 2961003



OUR TIP

Goggles in the cleanroom

Did you know? Cleanroom regulations recommend wearing goggles from class ISO 5 (GMP A/B) and higher.

The certified goggles from the PurGuard™ series are specifically developed for the use in cleanrooms and offer complete protection — of the product as well as person.

Soft frame materials guarantee the best wearer comfort even when worn for longer periods of time, while the latex-free headbands can be set individually and adjusted to the person using them.

Another highlight is the clever ventilation system of all PurGuard™ goggles. Large ventilation holes on the top and indirect vents on the bottom of the frame allow for a good airflow and reliably prevent the escaping of particles from the eye area.



CLEANROOM GOGGLES